


Ordering Information

Regenerative And Site Preparation



OSSEOTITE NT[®] Straight And Angled Osteotomes



NTOST0	Straight NT	Osteotome Kit
NTOST3	3.25mm Straight NT	Osteotome
NTOST4	4.0mm Straight NT	Osteotome
NTOST5	5.0mm Straight NT	Osteotome
NTOST6	6.0mm Straight NT	Osteotome
NTOST0A	Angled NT	Osteotome Kit
NTOST3A	3.25mm Angled NT	Osteotome
NTOST4A	4.0mm Angled NT	Osteotome
NTOST5A	5.0mm Angled NT	Osteotome
NTOST6A	6.0mm Angled NT	Osteotome

SUMMERS OSTEOTOMES



OST00	Summers Osteotome Kit	Includes #1-4
OST10	Summers Osteotome Kit	Includes #1-5
OST20	Summers Osteotome Kit	Includes #1-5 and FS
OST01	Summers Osteotome #1	Tip measures 1.6mm and flairs to 2.4mm at the 10mm depth mark. Prepares for the Osteotome #2
OST02	Summers Osteotome #2	Tip measures 1.9mm and flairs to 3.1mm at the 10mm depth mark. Prepares for the Osteotome #3.
OST03	Summers Osteotome #3	Tip measures 2.9mm and flairs to 3.4mm at the 10mm depth mark. Allows placement of a 3.75 or 4mm implant.
OST04	Summers Osteotome #4	Tip measures 3.2mm and flairs to 3.9mm at the 10mm depth mark. Allows placement of a 3.75 or 4mm implant.
OST05	Summers Osteotome #5	Tip measures 3.7mm and flairs to 4.9mm at the 10mm depth mark. Allows placement of a 5mm implant.
OSTFS	Summers Osteotome For Future Site	Tip measures 6mm and used when there is a low sinus floor.

SINUS ELEVATION



*SEKT1 Includes:

ME100	TIPL1
ME200	TIBC1
ME300	SEKTR
TIPS1	

SEKT1*	Sinus Elevation Kit
TIBC1	Titanium Nitride Coated Bone Carrier
ME100	Small Membrane Elevator
ME200	Medium Membrane Elevator
ME300	Large Membrane Elevator
TIPS1	Small Titanium Plugger
TIPL1	Large Titanium Plugger
SEKTR	Sinus Elevation Tray Only
TRE02	Trephine Bur 2mm I.D. 3mm O.D.
TRE04	Trephine Bur 4mm I.D. 5mm O.D.
TRE05	Trephine Bur 5mm I.D. 6mm O.D.
TRE06	Trephine Bur 6mm I.D. 7mm O.D.
TRE08	Trephine Bur 8mm I.D. 9.3mm O.D.

OSSIX™


6-Month Resorbable Collagen Membrane



OSX2530	25mm x 30mm GBR Membrane With Template 6-Month Resorbable Collagen Membrane
OSX3040	30mm x 40mm Large GBR Membrane With Template 6-Month Resorbable Collagen Membrane

BIOGRAN®

100% Synthetic Resorbable Bone Graft Material



SYRINGES*		
2100 - 0500	Syringe & 500mg Cup	1 Pack
2100 - 0750	Syringe & 750mg Cup	1 Pack
MIXING CUPS		
2100 - 0003	500mg Cups	7 Pack
2100 - 0004	750mg Cups	7 Pack
2100 - 0005	1500mg Cups	7 Pack

*Syringes are not pre-filled.

CALCIGEN™ ORAL

100% Synthetic Resorbable Bone Graft Binder/Stabilizer



CAL150

Calcigen Oral Bone Graft System

1.5 Gram Calcium Sulfate
.8 Gram Setting Solution

SAFESCRAPER®*

Disposable Cortical Bone Collector 3-Pack



3987

Curved Tube

3 Pack

3598

Straight Tube

3 Pack

*Not available in all regions.

PCCS® II

Platelet Concentrate Collection System



7426

PCCS II Centrifuge 120 Volt 50-60HZ

7427

PCCS II Centrifuge 240 Volt 50-60HZ

800-0550

PCCS II Disposable Kit (Sold in US only)

Contains The Following: One 60ml Syringe
One 10ml Syringe
One 18 Gauge Apheresis Needle Set
One 16 Gauge Angio Catheter
One Disposable Separation Tube
One Volume Gauge
One 50ml Bottle Of ACD-A

800-0570

PCCS II Disposable Kit (Sold outside US only)

Contains all items from kit above except for ACD-A

800-0575

PCCS II Disposable Kit (Sold outside US only)

Contains all items from Kit (800-0550) except with one 30ml bottle of ACD-A instead of a 50ml bottle

BONE AND SOFT TISSUE ALLOGRAFTS

For Periodontal And Oral Maxillofacial
Surgical Procedures



University of Miami Allograft* (DFDBA & FDBA)

0344	Crushed Cancellous	1/2gm	200-300µ
1449	Crushed Cancellous	1/2gm	300-500µ
1487**	Crushed Cancellous	4.0gm	300-500µ
0411	Crushed Cancellous	1/2gm	500-800µ
1354**	Crushed Cancellous	1.0gm	500-800µ
0605	Crushed Cortical	1/2gm	200-300µ
1450	Crushed Cortical	1/2gm	300-500µ
0393	Crushed Cortical	1/2gm	500-800µ
1451	Deminerlized Cancellous	1/2gm	200-300µ
1452	Deminerlized Cancellous	1/2gm	300-500µ
0400	Deminerlized Cancellous	1/2gm	500-800µ
0338	Deminerlized Cortical	1/2gm	200-300µ
1453	Deminerlized Cortical	1/2gm	300-500µ
0394	Deminerlized Cortical	1/2gm	500-800µ
1259**	Deminerlized Cortical	20cc	500-1000µ
0078	Dental Skin	10cm ²	
0079	Dental Skin	20cm ²	
0655	Dura Mater-Patch	30 cm ²	
0631	E.T.O. LT Hemi Mandible		
0632	E.T.O. Mandible		
0630	E.T.O. RT Hemi Mandible		
0716	Fascia Lata Patch	30cm ²	
0056	Illiic Crest Strip	Large	
0473	Illiic Crest Strip	Medium	
0055	Illiic Crest Strip	Small	
1454	Cancellous Powder	10cc	100-200µ
1455	Cancellous Powder	20cc	100-200µ
1456	Cancellous Powder	40cc	100-200µ
1457	Cancellous Powder	80cc	100-200µ
1458	Cortical Powder	10cc	100-200µ
1459	Cortical Powder	20cc	100-200µ
1460	Cortical Powder	40cc	100-200µ
1461	Cortical Powder	80cc	100-200µ
1462	Deminerlized Cancellous Powder	10cc	100-200µ
1463	Deminerlized Cancellous Powder	20cc	100-200µ
1464	Deminerlized Cancellous Powder	40cc	100-200µ
1465	Deminerlized Cortical Powder	10cc	100-200µ
1466	Deminerlized Cortical Powder	20cc	100-200µ
1467	Deminerlized Cortical Powder	40cc	100-200µ
1468	Deminerlized Cortical Powder	80cc	100-200µ



*Available for sale only in the United States.
Plus \$10.00 U of M handling charge.

**Special order, allow 4-6 weeks.

R. QUETIN BONE-MILL

Stainless Steel Bone Grinder



QBM001

Bone-Mill

Stainless Steel Bone Grinder

Implant Innovations, Inc.
Global Headquarters
4555 Riverside Drive
Palm Beach Gardens, FL 33410
1-800-342-5454
Outside The U.S.: +1-561-776-6700
Fax: +1-561-776-1272
www.3i-online.com



3i and design, Biogran, Implant Innovations, OSSEOTITE NT and PCCS are registered trademarks of Implant Innovations, Inc. Calcigen is a trademark of Biomet, Inc. OSSIX is a trademark of Colbar LifeScience Ltd. Safescraper is a registered trademark of CGM S.p.A Divisione Medicale META.
© 2005 Implant Innovations, Inc. All rights reserved.



ART745

REV F 08/05

SUBSIDIARIES

AUSTRALIA
Phone: +61-2-9855-4444
Fax: +61-2-9888-9900

BELGIUM
Phone: +32-2-5410290
Fax: +32-2-5410291

BRAZIL
Phone: +55-11-5081-4405
Fax: +55-11-5081-7484

CANADA
Phone: +514-956-9843
Fax: +514-956-9844

FRANCE
Phone: +33-1-41054343
Fax: +33-1-41054340

GERMANY
Phone: +49-721-255-1770
Fax: +49-721-255-1770

IRELAND
Phone: +35-31-477-3925
Fax: +35-31-402-9590

MEXICO
Phone: +52-55-5679-1619
Fax: +52-55-5684-8098

THE NETHERLANDS
Phone: +31-35-646-2830
Fax: +31-35-646-2666

NORDIC REGION
Phone: +46-40-17-6090
Fax: +46-40-17-6099

PORTUGAL
Phone: +351-21-000-1645
Fax: +351-21-000-1675

SPAIN
Phone: +34-93-470-59-50
Fax: +34-93-372-11-25

SWITZERLAND
Phone: +41-1-3804646
Fax: +41-1-3834655

U.K.
Phone: +44-1628-829314
Fax: +44-1628-820182

DISTRIBUTORS

ARGENTINA
Dentalmax, SA
Phone: +541-1482-71001
Fax: +541-1482-67373

AUSTRIA
Wieladent
Phone: +43-7672-93901
Fax: +43-7672-93903

CHILE
Cybel, SA
Phone: +56-2-2321883
Fax: +56-2-2330176

COLOMBIA
Implantes y Componentes
Phone: +571-612-9362
Fax: +571-620-5450

COSTA RICA
Implantec S.A.
Phone: +506-2-256411
Fax: +506-2-247620

EL SALVADOR
Dentimerc SA de CV
Phone: +503-263-6350
Fax: +503-263-6676

GREECE
Kostas Kornisorlis and Co.
Phone: +302310-269-079
Fax: +302310-555-573

ISRAEL
H.A. Systems
Phone: +972-3-6138777
Fax: +972-3-6138778

ITALY
Biomax, srl.
Phone: +39-0444-913410
Fax: +39-0444-913695

JAPAN
Implant Innovations Japan
Phone: +81-66-868-3012
Fax: +81-66-868-2444

KOREA
Jungsan Biomed Corp.
Phone: +82-2-516-1808
Fax: +82-2-514-9434

LEBANON
Tamer Freres s.a.l.
Phone: +961-1-485690
Fax: +961-1-510233

PANAMA
Odontomedica, S.A.
Phone: +507-2-239622
Fax: +507-2-239621

PARAGUAY
Andres H. Arce y Cia SRL
Phone: +595-21-208185
Fax: +595-21-496291

POLAND
Dental Depot
Phone: +48-71-341-3091
Fax: +48-71-343-6560

RUSSIA
Com-dental
Phone: +7095-797-6686
Fax: +7095-242-9567

SINGAPORE
Asia Implant Support & Services
Phone: +65-6223-2229
Fax: +65-6220-3538

TAIWAN
Kuo Hwa Dental Suppliers Co., Ltd.
Phone: +886-2-2226-1770
Fax: +886-2-2226-8747

THAILAND
PT Endeavour Co., Ltd.
Phone: +662-264-2574
Fax: +662-264-2573

UKRAINE
Olimpex, Ltd.
Phone: +38-044-209-5713
Fax: +38-044-482-3678

URUGUAY
Pro3Implant S.R.L.
Phone: +598-2-4034163
Fax: +598-2-4034163